

New replacement tubes from DUNLEE

AGITO Medical offers an extensive selection of new and refurbished medical imaging equipment and spare parts. In cooperation with DUNLEE, we are now able to offer you high quality CT replacement tubes with full OEM warranty.

Call +45 9635 0170 today or visit dunlee.com to learn more.

TubeType	Ordering No.	Focal Spot	kHU	Standard Warranty	
General Electric					
REEVO Series					
	REEVO240G This 8.0 MHU CT tube is designed as a replacement for use in GE Discovery™, GE Lightspeed™ and GE Optima™*	9896-055-96401	0.7 x 0.6 0.9 x 0.9	8,000	12 months, full
DA200 Series					
	DA200P40 This 6.3 MHU CT tube is designed as a replacement for use in GE Optima™*	4598-006-53451	0.9 x 0.7 1.2 x 1.1	6,300	12 months or 6,000 patient exams, full
	DA200ULTRA This 6.3 MHU CT tube is designed as a replacement for use in GE BrightSpeed™, GE Discovery™, GE Lightspeed™ and GE Optima™*	9896-055-88704	0.9 x 0.7 1.2 x 1.1	6,300	12 months, full
DA165 Series					
	DA165NP This 3.5 MHU CT tube is designed as a replacement for use in GE HiSpeed, GE BrightSpeed™ and GE Optima™*	9896-055-83901	0.9 x 0.7 1.2 x 1.2	3,500	14,000 slices/90 days full, 140,000 slices or 12 months pro-rated
DA135 Series					
	DA135 CT/e This 2.0 MHU CT tube is designed as a replacement for use in GE Brivo™, GE HiSpeed and GE CT Dual™*	9896-055-80800	0.7 x 0.6	2,000	8,000 slices/90 days full, 80,000 slices or 12 months pro-rated
Siemens					
S532/Akron Series					
	S532B/Akron B This 5.3 MHU CT tube is designed as a replacement for use in Siemens SOMATOM*	9896-055-95702	0.5 x 0.7 0.6 x 1.3	5,300	100%: 40,000 ss or 6 months 70%: 80,000 ss or 8 months 40%: 150,000 ss or 12 months
	S532Q/Akron Q This 5.3 MHU CT tube is designed as a replacement for use in Siemens SOMATOM*	9896-055-95902	0.5 x 0.7 0.8 x 1.2	5,300	100%: 30,000 ss or 6 months 70%: 60,000 ss or 8 months 40%: 120,000 ss or 12 months