

# System specifications

## CT10 – Toshiba Aquilion CXL – Trailer

<b>SYSTEM INFORMATION</b>	
Manufacturer and model	Toshiba Aquilion CXL
Number of slices	128 slices
Air or water cooling	Air-cooled
YOM	2005
<b>GANTRY AND DETECTOR</b>	
Gantry aperture	72 cm.
Detector model	Quantum Detector
<b>TUBE INFORMATION</b>	
Tube model	Tube Megacool 0.35s
<b>PATIENT TABLE</b>	
Max table load	205 kg.
Table length	188 cm.
Head rest incl.?	Yes
Positioning pads incl.?	Yes
<b>CONSOLE AND WORKSTATION(S)</b>	
Console model	
Console software version	V4:62GR017
Console software options	
Workstation (1) model	Vitrea Workstation
Workstation (1) SW version	
<b>INJECTOR SYSTEM</b>	
Manufacturer and model	Medrad Stellant Dual
Model no.	
<b>ACCESSORIES:</b>	