

# AGITO

M E D I C A L

Specification Sheet

## SIEMENS SOMATOM go.TOP

CT24 - Relocatable CT Solution



# SIEMENS SOMATOM go.TOP – Relocatable

This Relocatable Rental Solution from AGITO has SIEMENS SOMATOM go.Top CT equipment installed, which will provide the best scan for every type of patient – no matter the clinical demands and challenges. The CT scanner features a unique tablet-based mobile workflow, user guidance with SIEMENS GO technologies, and exclusive innovations such as Tin Filter low-dose technology.

This rental solution will meet the diverse needs of the modern imaging department with advanced clinical applications (such as advanced Cardio Package), with low patient dose and high throughput. Included are all the required accessories, with staff and patient comfort in mind, all housed inside a very high specification relocatable.



## System details

System ID	CT24	Dose Monitoring & Notific.	CARE KV, CARE Child, Dose Map
Type	Relocatable	Image Quality	Adaptive Signal Boost
Manufacturer & Model	SIEMENS Healthineers go.TOP	Metal Artifact Reduction	Check&Go Metal Detection, iMAR
No. of slices	128SL	Patient Monitoring	Patient CCTV installed
YOM (Unit)	2024	Cardiac Acquisition & Reconstruction	Advanced Cardiac Package (including Adaptive Cardio Sequence), Cardio Quick Sequence
Software Revision	VA30	ECG Trigger Monitor	PMM Included – Scan&Go Guidance on Electrode Positioning
Smart Tools & Accessories	Tablet & Remote Control	Max Patient Weight	227Kg
Gantry Aperture	70cm	Scannable Length	1,600 mm
Maximum Scan FOV	50cm	Table Feed Speed	1-200mm/s
Detector	STELLAR Detector	Vertical Travel Range	500mm - 885mm
Minimum Rotation Time	0.33s	Make / Model	BRACCO CTXpress
Tube	7.0MHU - X-Ray Tube	VPN / Firewall	OpenRad CUBE / Elite Networks
Generator	72KV	HIS/RIS & PACS	Easy Integration for full DICOM connectivity (via OpenRad CUBE)
KV Settings	70 - 140KV on 10KV steps	Remote Support	SRS Connectivity with OEM Service Contract & Remote Support
MA Settings	825 Ma (max)	Accessories	Coronal Supine Head Holder with positioning pads & bands, Table extension, Drip Stand, Body Mat & Bands, Knee Support wedge, Mattress Protector, Table Accessories set
SAFIRE	Sinogram Affirmed Iterative Reconstruction (CARE)		
Dose Reduction	ADMIRE (Advanced Modulated Iterative Reconstruction), FAST Planning, CARE DOSE 4D, X CARE, Flex DOSE Profile		

# SIEMENS SOMATOM go.TOP – Relocatable

Do you want to see how it looks inside the relocatable CT unit?

Take a 3D tour by scanning the QR code below.



## Unit features

Air Handling Unit - 10 Air Changes / Hour - Compliant to requirement for 'Interventional' procedures	Music System for Scan Room & Control Room
Electrically assisted door to scan room	Unit Hydraulic levelling System
BRACCO CT Exprés Dual Syringe Injector	Switchable Privacy Glass to Scan Room
Remote OEM service connectivity & OEM service contract on all equipment	Patient Changing Room with lockable door, cannulation chair, storage lockers and Patient Alarm
Sink with Hot & Cold Running Water in Scan Room	Ceiling mounted Patient Comfort Skylight in Scan room
Tri-zone air conditioning	Fire Alarm System – Fully Linkable to Hospital or Facility Alarm

## Unit details and dimensions

Type	Mobile
Length	14m
Width of unit	4m
Height	3.4m
Weight	36,600Kg (36T)
Min. area for placement	18m x 6m - Subject to Site Survey

## Electrical requirements from customer

400V, 3 phase, 50HZ, 200A. Marachell DS2 Female Connector with Dedicated Earth.
---

## Maintenance and service

System will be delivered with an OEM service and maintenance contract.

## Contact us

Mail: [rental@agitomedical.com](mailto:rental@agitomedical.com)

Phone: +45 96 35 01 70