

Philips Azurion 7 F12

Do you want to see how it looks inside the Relocatable Cath Lab unit?

Take a 3D tour by scanning the QR code below.



Unit features

FlexVision Pro – ability to view all connected applications at tableside	✓
StentBoost live - live enhanced visualization of device positioning and deployment	✓
FlexSpot reduce clutter and the number of separate workstations in control room	✓
Multimodality touch screen module ability to control all connected applications at tableside	✓
Dynamic coronary roadmap (unique to Philips) real time, automatic, motion compensated coronary imaging for guidance	✓

Unit features

Full Air Handling with heat exchange (25 air changes, H14 filters, +ve pressure)	✓
Medrad Mk7 injector	✓
Full network connectivity throughout	✓
Uninterruptible Power Supply (x-ray) (guaranteed 5 min fluoro)	✓
Isolated Power System with UPS (blue sockets throughout exam room)	✓
Unit fitted with hydraulic legs (Avoids requirement for crane – subject to site)	✓

System details

System ID	IT02-07
System model	Philips Azurion 7 F12
System YOM	2020/21
Max table load	250 kg
Single tube system	✓
HEMO with X3 monitoring	✓

Unit details and dimensions

Type	Modular
Length	15,00 m
Width	4,00 m
Height	3,30 m
Weight	28.000 kg
Min. area for placement	18,00 x 7,00 m

Maintenance and service

System will be delivered with a OEM service and maintenance contract.

Contact us

Mail: rental@agitomedical.com

Phone: +45 96 35 01 70