

# Philips Ingenia 1.5T

Powered by the breakthrough dStream architecture, the Ingenia 1.5T is the first MR system that brings MR signal digitization where it has never been before – in the RF coil, as close to the patient as possible. dStream unleashes the power of digitization by delivering a high purity MR signal for increased SNR, combined with Imaging Space enhanced workflow and ease of use for greater efficiency in your daily operations.



## System details

System ID	MR39
System model	Philips Ingenia 1.5T
System YOM	2014
No. of channels	Digital
Software rel.	5.3.1.2
Gantry opening	70 cm
Max table load	250 kg
Coils	Anterior, Posterior, 2 x Flex Small, NVC Head/Neck, NVC Head, Base, Knee 8ch., Foot/Ankle 8ch., Wrist 8ch., Sense Shoulder 8ch., Breast 7ch., Interfase

## Unit details and dimensions

Type	Trailer
Length	13,60 m
Width	2,55 m
Height	4,00 m
Weight	29.500 kg
Min. area for placement	16,60 x 7,60 m

## Unit features

Air-condition units in all rooms	✓
Coil storage cabinets	✓
Hydraulic patient lift	✓
Sound system w. CD/FM tuner	✓
IT connection in operator room	✓
Injector	✓
Changing area in control room	✓

## Maintenance and service

System will be delivered with a service and maintenance contract.

## Contact us

Mail: [rental@agitomedical.com](mailto:rental@agitomedical.com)

Phone: +45 96 35 01 70