

Siemens Aera 1.5T

Do you want to see how it looks inside the Mobile MRI unit?

Take a 3D tour by scanning the QR code below.



System details

System ID	MR40
System model	Siemens Aera 1.5T TIM 4G+DOT
System YOM	2021
No. of channels	24
Software rel.	Syngo XA30A
Gantry opening	70 cm
Max table load	250 kg
Coils	Head/Neck 16ch., Spine 24ch., Body 6ch., Flex Large 4ch., Flex Small 4ch., Flex Interface, Body 6ch. Long, Shoulder 16ch., Kit, Hand/Wrist 16ch., Tx/Rx Knee 15ch., Foot/Ankle 16ch..

Unit details and dimensions

Type	Trailer
Length	13,60 m
Width	2,55 m
Height	4,00 m
Weight	29.500 kg
Min. area for placement	16,60 x 7,60 m

Unit features

Air-condition units in all rooms	✓
Coil storage cabinets	✓
Hydraulic patient lift	✓
2x Sound system w. CD/DAP Radio for control & scan room	✓
IT connection in operator room	✓
Injector	✓
2xPatient changing area	✓

Maintenance and service

System will be delivered with a service and maintenance contract.

Contact us

Mail: rental@agitomedical.com

Phone: +45 96 35 01 70