

# Siemens Aera 1.5T

The Siemens Aera makes MRI exams easier and more comfortable than ever before. The roominess of the 70 cm Open Bore accommodates a large variety of patient sizes, shapes, and conditions. The ultra-short magnet (145 cm) allows many studies to be completed with the patient's head outside the bore.



## System details

System ID	MR28
System model	Siemens Aera 1.5T + DOT
System YOM	2019
No. of channels	48
Software rel.	SW Syngo MRE11C
Gantry opening	70 cm
Max table load	250 kg
Coils	Body long 18ch.,Shoulder Kit 16ch.,Hand/Wrist 16ch. Tx/Rx Knee Flare 15ch. Foot/Ankle 16ch.,Flex Large 4ch.,Flex Small 4ch.

## Unit details and dimensions

Type	Relocatable
Length	15,00 m
Width	4,20 m
Height	3,20 m
Weight	26.000 kg
Min. area for placement	17,00 x 6,20 m

## Unit features

Air-condition units in all rooms	✓
Coil storage cabinets	✓
Sound system w. CD/FM tuner	✓
IT connection in operator room	✓
Injector	✓
Patient changing area	✓

## Maintenance and service

System will be delivered with a service and maintenance contract.

## Contact us

Mail: [rental@agitomedical.com](mailto:rental@agitomedical.com)  
Phone: +45 96 35 01 70