

Canon Aquilion Prime SP

This relocatable Rental Solution from AGITO offers premium CT Scanning equipment which will be able to deliver even on the most challenging cases, from pediatric to bariatric & cardiac. The Aquilion Prime SP meets the diverse needs of the modern imaging department with advanced clinical applications, extremely low patient dose and high throughput.

Included are all the required accessories, with staff and patient comfort in mind, all housed inside a 'State of the Art' Relocatable.



System details

System ID	CT13/14
System model	Canon Aquilion Prime SP
System YOM	2021
No. of slices	160
Software rel.	10.4 SP008
Tube type	7.5MHU Megacool
Workstation	Vitrea Workstation
Gantry opening	78 cm
Max table load	220 kg

Unit details and dimensions

Type	Relocatable
Length	13,60 m
Width	3,60 m
Height	3,00 m
Weight	29.000 kg
Min. area for placement	20,00m x 6,00 m

Unit features

Triple Zone Air Conditioning & Staff Fridge	✓
Sink with Hot & Cold Running Water	✓
Switchable Privacy Glass	✓
Medrad Injector & contrast warmer	✓
Ceiling mounted Patient Comfort Skylight	✓
VISBION Cube for enhanced Connectivity Integration	✓
Patient changing aera	✓

Maintenance and service

System will be delivered with a service and maintenance contract.

Contact us

Mail: rental@agitomedical.com

Phone: +45 96 35 01 70