

# Toshiba Aquilion CXL

The Toshiba Aquilion CT scanner is designed for patient's safety and satisfaction. It has technology that meets your needs now and in the future. Advanced ergonomics are standard on every configuration. This system is capable of scanning a wide range of patients - from pediatric to bariatric. It also provides leading automated dose reduction capabilities to improve patient safety and increase your quality of care.



## System details

System ID	CT10
System model	Toshiba Aquilion CXL
System YOM	2005
No. of slices	128
Software rel.	V4:62GR017
Tube type	Tube Megacool 0.35s
Workstation	
Gantry opening	70 cm
Max table load	205 kg

## Unit details and dimensions

Type	Trailer
Length	13,50 m
Width	4,00 m
Height	4,00 m
Weight	27.000 kg
Min. area for placement	14,70 x 5,60 m

## Unit features

Air condition in all rooms	✓
Sink with running water	✓
IT connection in operator room	✓
Injector	✓
Sound system	✓
Coil/phantom storage cabinet	✓
Patient changing area	✓
Hydraulic lift	✓

## Maintenance and service

System will be delivered with a service and maintenance contract.

## Contact us

Mail: [rental@agitomedical.com](mailto:rental@agitomedical.com)

Phone: +45 96 35 01 70