

# GE Revolution EVO

Do you want to see how it looks inside the Mobile CT unit?

Take a 3D tour by scanning the QR code below.



## System details

System ID	CT17
System model	GE Revolution EVO 128
System YOM	2022
No. of slices	128
Software rel.	20HW16.10
Tube type	Performix 40 Plus
Minimum Rotation Time	0.35 s
Gantry opening	70 cm
Max table load	227 kg

## Unit details and dimensions

Type	Trailer
Length	13,56 m
Width	4,70 m
Height	4,00 m
Weight	27.000 kg
Min. area for placement	16,00 x 6,00 m

## Unit features

MEDRAD CENTARGO Injector (inbuilt contrast warmer)	✓
Advanced Cardio Package	✓
VISBION Cube for enhanced Connectivity Integration	✓
Patient Access Lift with Canopy & External Lighting	✓
Ceiling mounted Patient Comfort Skylight	✓
Switchable Privacy Glass	✓
Patient changing area	✓

## Maintenance and service

System will be delivered with a service and maintenance contract.

## Contact us

Mail: [rental@agitomedical.com](mailto:rental@agitomedical.com)

Phone: +45 96 35 01 70