

Siemens Aera 1.5T

The Siemens Aera makes MRI exams easier and more comfortable than ever before. The roominess of the 70 cm Open Bore accommodates a large variety of patient sizes, shapes, and conditions. The ultra-short magnet (145 cm) allows many studies to be completed with the patient's head outside the bore.



System details

System ID	MR34
System model	Siemens Aera 1.5T + DOT
System YOM	2013
No. of channels	48
Software rel.	VE11C
Gantry opening	70 cm
Max table load	250 kg
Coils	2 x body 18ch., Shoulder Large + Small, Wrist, TX/RX Knee, Foot/Ankle, PA36 Coil, Spine 32ch., Head/Neck 20Ch., Flex Small 4ch., Flex Large 4ch., Breast 18ch., TX/RX CP Head

Unit details and dimensions

Type	Relocatable
Length	15,00 m
Width	4,20 m
Height	3,30 m
Weight	31.000 kg
Min. area for placement	18,00 x 7,20 m

Unit features

Air-condition units in all rooms	✓
Coil storage cabinets	✓
Sound system w. CD/DAP Radio for control & scan room	✓
IT connection in operator room	✓
Injector	✓
Patient changing area	✓

Maintenance and service

System will be delivered with a service and maintenance contract.

Contact us

Mail: rental@agitomedical.com
Phone: +45 96 35 01 70