

Philips Ingenia 1.5T

Powered by the breakthrough dStream architecture, the Ingenia 1.5T is the first MR system that brings MR signal digitization where it has never been before – in the RF coil, as close to the patient as possible. dStream unleashes the power of digitization by delivering a high purity MR signal for increased SNR, combined with Imaging Space enhanced workflow and ease of use for greater efficiency in your daily operations.



System details

System ID	MR15
System model	Philips Ingenia 1.5T
System YOM	2015
No. of slices	Digital
Software rel.	Rel. 5.x
Gantry opening	70 cm
Max table load	250 kg
Coils	dStream Head/Neck Spine, Whole- Body, Flex S 2ch., Knee 16ch., Foot/ Ankle 8ch., Extremity S 8ch., Shoulder 8ch., Breast 7ch., Wrist 8ch., Microscopy Set

Unit details and dimensions

Type	Trailer
Length	13,60 m
Width	2,55 m
Height	4,00 m
Weight	29.500 kg
Min. area for placement	16,60 x 7,60 m

Unit features

Air-condition units in all rooms	✓
Coil storage cabinets	✓
Hydraulic patient lift	✓
Sound system w. CD/FM tuner	✓
IT connection in operator room	✓
Injector & Omega HP gradients	✓
Patient changing area	✓

Maintenance and service

System will be delivered with a service and maintenance contract.

Contact us

Mail: rental@agitomedical.com
Phone: +45 96 35 01 70