

Philips Incisive Premium

This Mobile CT Rental Solution from AGITO provides an operator work environment that is both flexible and easy to use, helping you to confidently achieve improved clinical outcomes, and enhance patient & staff experience.

The Philips Incisive Premium performs for you by allowing for intelligent decisions at every step. It provides applications that allow clinicians to improve workflow, and planning, as well as post processing analysis and review, at the same time ensuring extremely low patient dose. Included are all the required accessories, all housed inside a high-quality Mobile, that can be easily transported to your site.



System details

System ID	CT15
System model	Siemens Incisive Premium
System YOM	2021
No. of slices	128
Software rel.	V 5.0.0
Tube type	vMRC
Minimum Rotation Time	0.35 s
Gantry opening	72 cm
Max table load	205 kg

Unit details and dimensions

Type	Trailer
Length	13,56 m
Width	4,70 m
Height	4,00 m
Weight	27.000 kg
Min. area for placement	16,00 x 6,00 m

Unit features

MEDRAD CENTARGO Injector (inbuilt contrast warmer)	✓
Switchable Privacy Glass	✓
Ceiling mounted Patient Comfort Skylight	✓
Dual Zone Air Conditioning	✓
Music System in Scan, Control & Changing Room	✓
VISBION Cube for enhanced Connectivity Integration	✓
Patient changing area	✓

Maintenance and service

System will be delivered with a service and maintenance contract.

Contact us

Mail: rental@agitomedical.com
Phone: +45 96 35 01 70