

Siemens Avanto TIM 1.5T

The Siemens Magnetom Avanto 1.5 Tesla MR machine is a popular, affordable MRI scanner. Its short bore allows it to fit in tight rooms or where a compact MRI is a more effective solution. The system features syngo SPACE, which replaces time-consuming multiple 2D acquisitions with a single 3D image, allowing you to visualize very small details.



System details

System ID	MR25
System model	Siemens Avanto 1.5T
System YOM	2009
No. of channels	18
Software rel.	MR B19
Gantry opening	60 cm
Max table load	250 kg
Coils	Flex Large + Small 4ch., Interface 4ch., High Re-solution Knee Array 8ch, Body Matrix, Coil Extremity, Shoulder Array, Head Matrix, Neck Matrix, Spine Matrix

Unit details and dimensions

Type	Trailer
Length	13,55 m
Width	2,55 m
Height	4,00 m
Weight	28.000 kg
Min. area for placement	16,55 x 5,55 m

Unit features

Air-condition units in all rooms	✓
Coil storage cabinets	✓
Sound system w. CD/FM tuner	✓
IT connection in operator room	✓
Injector	✓

Maintenance and service

System will be delivered with a service and maintenance contract.

Contact us

Mail: rental@agitomedical.com
Phone: +45 96 35 01 70