

Siemens Aera 1.5T

The Siemens Aera makes MRI exams easier and more comfortable than ever before. The roominess of the 70 cm Open Bore accommodates a large variety of patient sizes, shapes, and conditions. The ultra-short magnet (145 cm) allows many studies to be completed with the patient's head outside the bore.



System details

System ID	MR26
System model	Siemens Aera 1.5T + DOT
System YOM	2019
No. of channels	24
Software rel.	Syngo MRE11C
Gantry opening	70 cm
Max table load	250 kg
Coils	Flex Large 4ch., Flex Small 4ch., Body 6ch. long Shoulder Large+Small 16ch., Hand/Wrist 16ch. Tx/Rx Knee 15ch. Flare. Foot/Ankle 16ch., Spine 24ch. Head/Neck 16ch..

Unit details and dimensions

Type	Trailer
Length	13.60 m
Width	2.55 m
Height	4.00 m
Weight	29.500 kg
Min. area for placement	16,60 x 7,60 m

Unit features

Air-condition units in all rooms	✓
Coil storage cabinets	✓
Hydraulic patient lift	✓
Sound system w. CD/FM tuner	✓
IT connection in operator room	✓
Injector	✓
Patient changing area	✓

Maintenance and service

System will be delivered with a service and maintenance contract.

Contact us

Mail: rental@agitomedical.com
Phone: +45 96 35 01 70