

# GE Revolution EVO

The Revolution EVO has a continuously rotating x-ray source and detector array. It can be used as a whole-body CT scanner for general imaging use, and it offers options for cardiac, pediatric, oncology, vascular, orthopedic, and neuro imaging.

The Smart Flow Technologies enable hands-free patient positioning, exam prescription from the patient's side, integrated injections, real-time reconstruction and access to advanced applications on the console.



## System details

System ID	CT09
System model	GE Revolution EVO
System YOM	2020
No. of slices	128
Software rel.	NIO64
Tube type	Performix 40 Plus
Workstation	Volumeshare 7
Gantry opening	70 cm
Max table load	306 kg

## Unit details and dimensions

Type	Relocatable
Length	13,00 m
Width	4,55 m
Height	3,20 m
Weight	29.000 kg
Min. area for placement	13,00 x 4,55 m

## Unit features

Air condition in all rooms	✓
Sink with running water	✓
IT connection in operator room	✓
Injector	✓
Sound system	✓

## Maintenance and service

System will be delivered with a service and maintenance contract.

## Contact us

Mail: [rental@agitomedical.com](mailto:rental@agitomedical.com)  
Phone: +45 96 35 01 70