

# Siemens Aera 1.5T

The Siemens Aera makes MRI exams easier and more comfortable than ever before. The roominess of the 70 cm Open Bore accommodates a large variety of patient sizes, shapes, and conditions. The ultra-short magnet (145 cm) allows many studies to be completed with the patient's head outside the bore.



## System details

System ID	MR43
System model	Siemens Aera 1.5T TIM 4G+DOT
System YOM	2021
No. of channels	24
Software rel.	Syngo XA30A
Gantry opening	70 cm
Max table load	250 kg
Coils	Head/Neck 16ch., Spine 24ch., Body 6ch., Flex Large 4ch., Flex Small 4ch., Flex Interface, Body 6ch. Long, Shoulder 16ch., Kit, Hand/Wrist 16ch., Tx/Rx Knee 15ch., Foot/Ankle 16ch..

## Unit details and dimensions

Type	Trailer
Length	13,60 m
Width	2,55 m
Height	4,00 m
Weight	29.500 kg
Min. area for placement	16,60 x 7,60 m

## Unit features

Air-condition units in all rooms	✓
Coil storage cabinets	✓
Hydraulic patient lift	✓
2x Sound system w. CD/DAP Radio for control & scan room	✓
IT connection in operator room	✓
Injector	✓
2x Patient changing area	✓

## Maintenance and service

System will be delivered with a service and maintenance contract.

## Contact us

Mail: [rental@agitomedical.com](mailto:rental@agitomedical.com)  
Phone: +45 96 35 01 70