

Philips Intera Achieva 1.5T

The Philips Achieva 1.5T helps you to do more in every way. Fast exams, a wide application range, and superb workflow are just some of the results of the Achieva's do more philosophy. This system provides exceptional quality in challenging applications, like large field-of-view abdominal imaging, uniform fat suppression in virtually all anatomies, and consistent off-center imaging, across all patients.



System details

System ID	MR24
System model	Philips Intera Achieva 1.5T
System YOM	2007 (upgr. 2011)
No. of channels	24
Software rel.	Rel. 3.2.3.4
Gantry opening	60 cm
Max table load	250 kg
Coils	Syn Spine, Syn Body, Quadrature Knee/Foot, QHC El Load ASC-NT, NV 16ch., Sense XL Torso 16 ch., Sense Flex Medium, Sense Knee 8ch.

Unit details and dimensions

Type	Trailer
Length	13,55 m
Width	2,55 m
Height	4,00 m
Weight	28.000 kg
Min. area for placement	16,55 x 5,55 m

Unit features

Air-condition units in all rooms	✓
Coil storage cabinets	✓
Hydraulic patient lift	✓
Sound system w. CD/FM tuner	✓
IT connection in operator room	✓
Injector	✓
DVD-RW writer	✓

Maintenance and service

System will be delivered with a service and maintenance contract.

Contact us

Mail: rental@agitomedical.com

Phone: +45 96 35 01 70